

TSKgel SP-STAT and TSKgel CM-STAT columns

Non-Porous Cation Exchange Columns for High Speed and High Resolution Analysis of Biomolecules

TSKgel
PRODUCT OVERVIEW


Introduction

TSKgel SP-STAT and TSKgel CM-STAT cation exchange columns allow fast equilibration and analysis, as well as isolation, of complex biomolecules. Both TSKgel columns are packed with 7 or 10µm mono-disperse, non-porous resin particles of which the surface consists of an open access network of multi-layered cation exchange groups (see Figure 1). The innovative bonding chemistry, combined with a relatively large particle size, result in a respectable loading capacity and a low operating pressure, attributes not found in traditional mono-disperse, non-porous resins.

Product Highlights

- Very efficient chromatography for high as well as low MW solutes made possible by novel bonding chemistry and the absence of micro-pores
- High speed and high resolution analysis of biomolecules
- Higher adsorption capacities and lower pressures compared with competitive non-porous columns
- 7 or 10µm particles for SP and CM chemistries

Figure 1.




Applications

Fast Separations

The fast separation of protein standards was investigated using short cation exchange columns (see Figure 2). A TSKgel SP-STAT column shows superior resolution, better peak shape, and a shorter analysis time (< 60 seconds) compared to a competitive monolithic SP-type column.

Figure 2.



Reaction Monitoring


A sample of β -lactoglobulin (5mg/mL) was reacted with polyethylene glycol (5kDa) in a pH 6.5 phosphate buffer. The formation of pegylated protein reaction products was monitored in 5 minute intervals on a 3.5cm TSKgel SP-STAT column. As demonstrated in Figure 3, peak areas of mono-, di-, and tri-pegylated β -lactoglobulin increased with reaction time, while the area of unreacted β -lactoglobulin declined.



TOSOH BIOSCIENCE

TOSOH


Figure 3.



Antibody Analysis

The analysis profiles for five antibodies separated on a TSKgel CM-STAT column were compared with the profiles obtained on a competitive WCX-10 column (Figure 4). Similar or higher resolution profiles were obtained on TSKgel CM-STAT in approximately half the time.

Figure 4.



Ordering Information

Part #	Description	Matrix	Housing	ID (mm)	Length (cm)
21963	TSKgel SP-STAT, 10 μ m	Polymer	Stainless Steel	3	3.5
21964	TSKgel SP-STAT, 7 μ m	Polymer	Stainless Steel	4.6	10
21965	TSKgel CM-STAT, 10 μ m	Polymer	Stainless Steel	3	3.5
21966	TSKgel CM-STAT, 7 μ m	Polymer	Stainless Steel	4.6	10



TOSOH BIOSCIENCE

TOSOH Bioscience LLC
 3604 Horizon Drive, Suite 100
 King of Prussia, PA 19406
 Orders & Service: (800) 366-4875
 Fax: (610) 272-3028
www.separations.us.tosohbioscience.com
 email: info.tbl@tosoh.com