



Characterization of Biomacromolecules by FPLC coupled with the LenS₃ MALS Detector

Introduction

Size exclusion chromatography (SEC), specifically gel filtration chromatography (GFC), is an important tool for the separation of biomaterials and biopharmaceutical products such as monoclonal antibodies (mAbs). With the increasing need to acquire an in-depth understanding of physical properties, the addition of a multi-angle light scattering (MALS) detector to the common chromatography system has attracted great attention.

ÄKTA pure[®] is a commonly utilized FPLC system for the purification of proteins, peptides, and nucleic acids. This application note demonstrates the compatibility of the ÄKTA pure with the LenS₃ MALS detector. The coupling provides the user with the capability not only to qualitatively detect each fraction that is separated by the column connected to ÄKTA system, but also to quantitatively determine the absolute molecular weight and the size of each species, including the targeted molecules, impurities, aggregates, and fractions.

Material and Methods

Instrumentation: ÄKTA pure 25 in-line with UV detector (@ 280 nm) and LenS₃ MALS detector
 Column: 1 × Superdex[®] 200 Increase 10/300 GL, 8.6 μm, 10 mm ID × 30 cm
 Mobile phase: 100 mmol/L sodium phosphate, 100 mmol/L sodium sulfate, 0.01% sodium azide, pH 6.8
 Flow rate: 0.60 mL/min
 Temperature: ambient
 Injection volume: 35 μL

Results and Discussion

Monoclonal antibodies (mAbs)

Samples & Conditions

- NIST mAb (0.48 mg/mL); Trastuzumab (0.79 mg/mL, Herceptin[®] biosimilar); Adalimumab (0.45 mg/mL, Humira[®] biosimilar)
- $dn/dc = 0.187$ mL/g; dA/dc (NIST mAb) = 1.42 mL/g; dA/dc (trastuzumab) = 1.48; dA/dc (adalimumab) = 1.4

Figure 1. Molecular weight distribution overlaid with UV signals of NIST mAb (red), trastuzumab (blue), and adalimumab (green).




Table 1. Molecular weight characterization of mAbs by LenS₃

Peak	Fragment	Monomer	Dimer	Trimer
NIST mAb	N/A	147,198	N/A	N/A
Trastuzumab	116,081 (1.40%)	147,118 (79.62%)	306,051 (12.40%)	487,714 (0.64%)
Adalimumab	97,930 (2.72%)	146,059 (96.53%)	329,747 (0.75%)	N/A

Conclusions


- The calculated molecular weights of fragments, monomer, dimer and trimer from different mAbs show that these molecular weights are very similar to the reported molecular weight from literature. Results demonstrated that LenS₃ can determine mAb molecular weights accurately.
- LenS₃ provides an accurate and absolute MW calculation that is independent of retention time shift, caused e.g. by column use history.

Bovine Serum Albumin (BSA)

Samples & Conditions

- BSA from Sigma Aldrich (purity 98% by agarose gel electrophoresis)
- Concentration = 4.68 mg/mL
- dn/dc = 0.185 mL/g; dA/dc = 0.667 mL/g

➤ **Figure 2.** LALS (red), RALS (blue), and UV (yellow) signals of the BSA, and its molecular weight distribution (green).



➤ **Table 2.** Molecular weight characterization of BSA by LenS₃

Peak	Assignment	Retention time (min)	MW (Da)	%
1	Monomer	22.890	66,753	74.39
2	Dimer	19.625	130,473	15.56
3	Trimer	17.786	197,473	4.94
4	Tetramer	16.611	263,962	3.53
5	Pentamer	16.056	338,954	1.14

➤ **Table 3.** Reproducibility of molecular weight characterization of BSA by LenS₃

Injection #	Monomer	Dimer	Trimer	Tetramer	Pentamer
1	66,753	130,473	197,473	263,962	338,954
2	66,902	130,207	196,442	261,635	339,866
3	66,979	130,927	196,226	265,314	341,306
Average	66,878	130,536	196,714	264,637	340,042
%RSD	0.172	0.279	0.339	0.706	0.349

Conclusions


- LenS₃ has great sensitivity to detect and calculate not only the protein monomer but also the aggregates even at trace amounts (e.g. pentamer with only 1.14%)!
- The sensitivity and robust performance of LenS₃ enable achieving a %RSD <1% even for the pentamer.

Detection Linearity

Samples & Conditions

- BSA from Sigma Aldrich (purity 98% by agarose gel electrophoresis) with various injection amounts.

➤ **Figure 3.** RALS signal overlay of BSA with various injection loads.



Conclusions

- LenS₃ provides an exceptional detection linearity to different sample loading amounts (below 0.5% deviation).

Summary

Coupled with an ÄKTA pure FPLC system, the LenS₃ MALS detector provides:

- Unprecedented sensitivity for characterization of protein monomer and aggregates.
- Superb reproducibility and precision.
- Accurate and reliable molecular weight calculation (retention time shift independent!).
- Exceptional linearity to different injection loading amounts.

This complimenting use of the LenS₃ MALS detector with an ÄKTA pure FPLC system has not been approved or sponsored by Cytiva. Tosoh Bioscience does not have a relationship with Cytiva.

TSKgel and Tosoh Bioscience are registered trademarks of Tosoh Corporation.
 LenS is a registered trademark of Tosoh Bioscience LLC in the USA, India and Japan.
 SECview is a registered trademark of Tosoh Bioscience LLC in the USA, EU, and India; and of Tosoh Corporation in Japan.
 ÄKTA and Superdex are the registered trademarks of Cytiva.
 Herceptin is a registered trademark of Genentech, Inc.
 Humira is a registered trademark of AbbVie Biotechnology Ltd.