

TSKgel® FcR Products

Part Number:	23513, TSKgel FcR-III-NPR (non-porous resin), 5 µm, 4.6 mm ID × 7.5 cm
	23532, TSKgel FcR-III-5PW, 10 µm, 7.8 mm ID × 7.5cm

This sheet contains the recommended operating conditions and the specifications for TSKgel FcR columns. Installation instructions and column care information are described in a separate Instruction Manual.

A. OPERATING CONDITIONS	
1. Shipping Solvent:	0.65 mmol/L Critic acid, 9.35 mmol/L Trisodium citrate containing 0.025% ProClin® 300, pH 6.5
2. Standard Flow Rate:	1.0 mL/min (4.6 mm ID) 0.3 mL/min (7.8 mm ID)
3. Max Flow Rate:	1.0 mL/min (4.6 mm ID) 1.2 mL/min (7.8 mm ID)
NOTE:	When a buffer with high viscosity is used, the maximum flow rate may have to be reduced so as not to exceed the maximum pressure drop. When changing solvents, use a flow rate equal to 25% of the maximum flow rate.
4. Max. Pressure:	9.0 MPa (4.6 mm ID) 1.0 MPa (7.8 mm ID)
5. Temperature:	15 – 25 °C.
6. pH Range:	4.0 – 8.0
7. Organic Concentration:	TSKgel FcR-III-NPR is a protein ligand immobilized resin so that any organic solvent and denaturing agent like SDS, urea or guanidine must not be used as mobile phase.
8. Cleaning Solvents:	1. Turn the column in reverse flow direction and run at half the maximum flow rate. 2. Make 3-5 injections of 0.5 - 2 mL of a buffer containing 0.5 mol/L NaCl or 20 % ethanol 3. Turn column in normal flow direction and equilibrate in mobile phase for at least 45 minutes
9. Storage:	1. Procedure: a. Replace the column contents with the shipping solvent, disconnect the column from the instrument, seal both ends with the end plugs, and store. b. After disconnecting the column from the instrument, wash the instrument tubing with distilled water or ion exchange water. Use the solvent replacement flow rate during cleaning and when replacing with the shipping solvent. 2. Storage temperature: 2 to 8 °C
NOTE:	

B. SPECIFICATIONS	
The performance of TSKgel FcR columns is tested under the conditions described in the data sheet. All columns have passed the following quality control specifications:	
1. Number of Theoretical Plates (N):	≥ 170 (4.6 mm ID) ≥ 2100 (7.8 mm ID)
2. Asymmetry (AF) / Symmetry (S) Factor:	AF: 1.0 - 1.8 (4.6 mm ID) S: 0.8 - 1.6 (7.8 mm ID)

DS1243 22MARCH2021

www.tosohbioscience.com

Tosoh Bioscience LLC
3604 Horizon Drive, Suite 100
King of Prussia, PA 19406
Phone: (484) 805-1219
Orders & Technical Service: (800) 366-4875
FAX: (610) 272-3028

Tosoh Bioscience GmbH
Im Leuschnerpark 4, 64347
Griesheim, Germany
Phone/ FAX: +49-6155-7043700