

# TOYOPEARL® Sulfate-650F

Strong cation exchange resin for capture and removal of mAb aggregates


TOYOPEARL  
PRODUCT OVERVIEW

## Introduction

TOYOPEARL Sulfate-650F resin is a strong cation exchange resin that exhibits high salt tolerance. This resin offers the strongest capture of monoclonal antibody (mAb) aggregates over a wide pH range without losing its binding capacity for mAb. With the use of optimized binding conditions, a dynamic binding capacity of >120 g/L of mAb can be easily achieved with TOYOPEARL Sulfate-650F. This high binding capacity translates into lower operating costs per gram of antibody produced.

A TOYOPEARL HW-65F polymeric bead has been functionalized with a sulfate (SO<sub>4</sub><sup>-</sup>) group (see *Structure*). The 100 nm pore size of this resin, along with proprietary bonding technology, makes TOYOPEARL Sulfate-650F ideal for applications performed in physiological conditions or for post-protein A removal of aggregates.

STRUCTURE:




## Product Attributes

Pore size (mean):	100 nm
Particle size (mean):	45 µm (F-grade)
Pressure rating:	0.3 MPa
Shipping buffer:	20% ethanol
pH stability:	2-13

## High Dynamic Binding Capacity

TOYOPEARL Sulfate-650F offers high dynamic binding capacities for IgG, obtainable even at higher flow rates, as shown in *Figure 1*.


**Figure 1.** Dynamic Binding Capacities



## Effectively Capture and Remove Aggregates from mAb Sample

TOYOPEARL Sulfate-650F resin is effective at removing aggregates from mAbs, as demonstrated in *Figure 2*. A protein A-purified IgG sample was loaded onto a TOYOPEARL Sulfate-650F column, fractions were collected using an ÄKTA® and further analyzed using a TSKgel® G3000SW<sub>XL</sub> HPLC column. TOYOPEARL Sulfate-650F resin provides stronger binding of mAb aggregates, resulting in the high resolution separation of monomer and aggregates.

**Figure 2.** Aggregate Removal




TOSOH BIOSCIENCE

TOSOH

Figure 3 demonstrates the excellent pressure-flow rate properties of the TOYOPEARL Sulfate-650F resin. A flow rate of >600 cm/hr on a large process column is easily achieved at a pressure drop of only 0.2 MPa.


Figure 3. Pressure-Flow Rate Curve on Large Process Column (30 cm ID x 20 cm bed height)



## Alkaline Stability

TOYOPEARL Sulfate-650F resin is stable in 0.5 mol/L NaOH (Figure 4). It can be stored in this solution for up to 8 weeks without loss in its binding capacity.

Figure 4. Alkaline Stability of TOYOPEARL Sulfate-650F



## Ordering Information

Part #	Description	Resin Volume
23467	TOYOPEARL Sulfate-650F	100 mL
23468	TOYOPEARL Sulfate-650F	250 mL
23469	TOYOPEARL Sulfate-650F	1 L
23470	TOYOPEARL Sulfate-650F	5 L
23471	TOYOPEARL Sulfate-650F	50 L
23472	TOYOPEARL Sulfate-650F	1 mL x 6
23473	TOYOPEARL Sulfate-650F	5 mL x 6
45117	TOYOPEARL Sulfate-650F MiniChrom Column	5 mL
45027	ToyoScreen Sulfate-650F RoboColumn	8 x 200 µL
45028	ToyoScreen Sulfate-650F RoboColumn	8 x 600 µL
OC41MDSLFT-650F	TOYOPEARL Sulfate-650F Resin Seeker 96-Well Plate	20 µL resin bed

TOYOPEARL, TSKgel and Tosoh Bioscience are registered trademarks of Tosoh Corporation  
ÅKTA is a registered trademark of GE Healthcare  
AxiChrom is a registered trademark of GE Healthcare Bioprocess R&D AB